

## Zinger Chair Bottle Holder Installation Instructions

---

1. Open the clamp by unthreading the lever until the jaws of the clamp fit over the silver front tube of your Zinger

---



If you have any questions please call us at 800-958-8324

---

2. Close the jaws only a little by turning the lever to the right about three full turns, then fold the lever down as shown below. It should require moderate force to fold it down

---



---

If the clamp is not tight enough, open the lever and rotate it one more turn to the right.

---

If the clamp is too tight to fold the lever down, rotate it one turn to the left, and then try folding it down again.

---

## Zinger Chair Bottle Holder Installation Instructions

---

1. Open the clamp by unthreading the lever until the jaws of the clamp fit over the silver front tube of your Zinger

---



If you have any questions please call us at 800-958-8324

---

2. Close the jaws only a little by turning the lever to the right about three full turns, then fold the lever down as shown below. It should require moderate force to fold it down

---



---

If the clamp is not tight enough, open the lever and rotate it one more turn to the right.

---

If the clamp is too tight to fold the lever down, rotate it one turn to the left, and then try folding it down again.

---